

Product Datasheet

Mouse PDCD1LG1 ELISA Kit (Ready to Use) (orb553956)

Description	Mouse Programmed Cell Death Protein 1 Ligand 1 (PDCD1LG1) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9EP73
Sensitivity	0.058ng/mL
Entrez	60533
Expiration Date	Please enquire.