

Product Datasheet

Mouse PROZ ELISA Kit (Ready to Use) (orb553988)

Description	Mouse Protein Z (PROZ) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability plasma.
Sample Types	plasma.
Uniprot ID	Q9CQW3
Sensitivity	6.5pg/mL
Entrez	66901
Expiration Date	Please enquire.