

Product Datasheet

Mouse PKM2 ELISA Kit (Ready to Use) (orb553997)

Description	Mouse Pyruvate Kinase, Muscle (PKM2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents ...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P52480
Sensitivity	0.144ng/mL
Entrez	18746
Expiration Date	Please enquire.