

Product Datasheet

Mouse RSP01 ELISA Kit (Ready to Use) (orb553999)

Description	Mouse R-Spondin 1 (RSP01) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z132
Sensitivity	5.8pg/mL
Entrez	192199
Expiration Date	Please enquire.