

Product Datasheet

Mouse REG3a ELISA Kit (Ready to Use) (orb554007)

Description	Mouse Regenerating Islet Derived Protein 3 Alpha (REG3a) Sandwich ELISA KitReady to Use: Pre-diluted...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	009037
Sensitivity	6.3pg/mL
Entrez	19694
Expiration Date	Please enquire.