

Product Datasheet

Mouse SOSTDC1 ELISA Kit (Ready to Use) (orb554041)

Description	Mouse Sclerostin Domain Containing Protein 1 (SOSTDC1) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9CQN4
Sensitivity	0.068ng/mL
Entrez	66042
Expiration Date	Please enquire.