

## Product Datasheet

### Mouse SEMA3E ELISA Kit (Ready to Use) (orb554057)

<b>Description</b>	Mouse Semaphorin 3E (SEMA3E) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conv...
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P70275</a>
<b>Sensitivity</b>	0.122ng/mL
<b>Entrez</b>	<a href="#">20349</a>
<b>Expiration Date</b>	Please enquire.