

Product Datasheet

Mouse SIRT6 ELISA Kit (Ready to Use) (orb554085)

Description	Mouse Sirtuin 6 (SIRT6) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenien...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P59941
Sensitivity	0.063ng/mL
Entrez	50721
Expiration Date	Please enquire.