

Product Datasheet

Mouse SMOX ELISA Kit (Ready to Use) (orb554102)

Description	Mouse Spermine Oxidase (SMOX) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q99K82
Sensitivity	0.122ng/mL
Entrez	228608
Expiration Date	Please enquire.