

Product Datasheet

Mouse SUOX ELISA Kit (Ready to Use) (orb554122)

Description	Mouse Sulfite Oxidase (SUOX) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8R086
Sensitivity	0.054ng/mL
Entrez	211389
Expiration Date	Please enquire.