

Product Datasheet

Mouse SPD ELISA Kit (Ready to Use) (orb554133)

Description	Mouse Surfactant Associated Protein D (SPD) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, lung lavage fluid, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates and other biological fluids.
Uniprot ID	P50404
Sensitivity	0.285ng/mL
Entrez	20390
Expiration Date	Please enquire.