

Product Datasheet

Mouse TLR2 ELISA Kit (Ready to Use) (orb554178)

Description	Mouse Toll Like Receptor 2 (TLR2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
Reactivity	Mouse
Range	3.12-200pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9QUN7
Sensitivity	1.42pg/mL
Entrez	24088
Expiration Date	Please enquire.