

## Product Datasheet

### Mouse TLR4 ELISA Kit (Ready to Use) (orb554180)

<b>Description</b>	Mouse Toll Like Receptor 4 (TLR4) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9QUK6</a>
<b>Sensitivity</b>	0.065ng/mL
<b>Entrez</b>	<a href="#">21898</a>
<b>Expiration Date</b>	Please enquire.