

Product Datasheet

Mouse TFR2 ELISA Kit (Ready to Use) (orb554188)

Description	Mouse Transferrin Receptor 2 (TFR2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	Mouse
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9JKX3
Sensitivity	0.61ng/mL
Entrez	50765
Expiration Date	Please enquire.