

Product Datasheet

Mouse TGFb2 ELISA Kit (Ready to Use) (orb554189)

Description	Mouse Transforming Growth Factor Beta 2 (TGFb2) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P27090
Sensitivity	13.2pg/mL
Entrez	21808
Expiration Date	Please enquire.