

Product Datasheet

Mouse TRPM4 ELISA Kit (Ready to Use) (orb554196)

Description	Mouse Transient Receptor Potential Cation Channel Subfamily M, Member 4 (TRPM4) Sandwich ELISA KitRe...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q7TN37
Sensitivity	0.048ng/mL
Entrez	68667
Expiration Date	Please enquire.