

Product Datasheet

Mouse TRPV2 ELISA Kit (Ready to Use) (orb554200)

Description	Mouse Transient Receptor Potential Cation Channel Subfamily V, Member 2 (TRPV2) Sandwich ELISA KitRe...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9WTR1
Sensitivity	0.049ng/mL
Entrez	22368
Expiration Date	Please enquire.