

Product Datasheet

Mouse TFF2 ELISA Kit (Ready to Use) (orb554206)

Description	Mouse Trefoil Factor 2 (TFF2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9QX97
Sensitivity	12.5pg/mL
Entrez	21785
Expiration Date	Please enquire.