

Product Datasheet

Mouse TNFRSF17 ELISA Kit (Ready to Use) (orb554229)

Description	Mouse Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17) Sandwich ELISA Kit Ready to Us...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O88472
Sensitivity	0.101ng/mL
Entrez	21935
Expiration Date	Please enquire.