

Product Datasheet

Mouse Tie1 ELISA Kit (Ready to Use) (orb554236)

Description	Mouse Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) Sandwich ELISA ...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q06806
Sensitivity	0.247ng/mL
Entrez	21846
Expiration Date	Please enquire.