

Product Datasheet

Mouse VIP ELISA Kit (Ready to Use) (orb554263)

Description	Mouse Vasoactive Intestinal Peptide (VIP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	7.81-500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P32648
Sensitivity	2.26pg/mL
Entrez	22353
Expiration Date	Please enquire.