

Product Datasheet

Mouse DBP ELISA Kit (Ready to Use) (orb554276)

Description	Mouse Vitamin D Binding Protein (DBP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
Reactivity	Mouse
Range	6.25-400ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P21614
Sensitivity	2.27ng/mL
Entrez	14473
Expiration Date	Please enquire.