

## Product Datasheet

### Mouse VWF ELISA Kit (Ready to Use) (orb554280)

<b>Description</b>	Mouse Von Willebrand Factor (VWF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
<b>Reactivity</b>	Mouse
<b>Range</b>	0.25-16ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability plasma.
<b>Sample Types</b>	plasma.
<b>Uniprot ID</b>	<a href="#">Q8CIZ8</a>
<b>Sensitivity</b>	0.11ng/mL
<b>Entrez</b>	<a href="#">22371</a>
<b>Expiration Date</b>	Please enquire.