

Product Datasheet

Porcine aHSG ELISA Kit (Ready to Use) (orb554301)

Description	Porcine Alpha-2-Heremans Schmid Glycoprotein (aHSG) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Porcine
Range	15.6-1000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P29700
Sensitivity	6.7ng/mL
Expiration Date	Please enquire.