

Product Datasheet

Porcine Angiopoietin 2 ELISA Kit (Ready to Use) (orb554305)

Description	Porcine Angiopoietin 2 (Angiopoietin 2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Porcine
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9BDY7
Sensitivity	0.54ng/mL
Entrez	396730
Expiration Date	Please enquire.