

## Product Datasheet

# Porcine APOA2 ELISA Kit (Ready to Use) (orb554309)

<b>Description</b>	Porcine Apolipoprotein A2 (APOA2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
<b>Reactivity</b>	Porcine
<b>Range</b>	4.68-300µg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<a href="#">F1S1A9</a>
<b>Sensitivity</b>	1.2µg/mL
<b>Entrez</b>	<a href="#">100153243</a>
<b>Expiration Date</b>	Please enquire.