

Product Datasheet

Porcine BNP ELISA Kit (Ready to Use) (orb554319)

Description	Porcine Brain Natriuretic Peptide (BNP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagen...
Reactivity	Porcine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P07634
Sensitivity	13.2pg/mL
Entrez	396844
Expiration Date	Please enquire.