

Product Datasheet

Porcine CNTF ELISA Kit (Ready to Use) (orb554327)

Description	Porcine Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection rea...
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	F1RMF2
Sensitivity	22.4pg/mL
Entrez	100518821
Expiration Date	Please enquire.