

Product Datasheet

Porcine C3 ELISA Kit (Ready to Use) (orb554330)

Description	Porcine Complement Component 3 (C3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	Porcine
Range	0.39-25ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01025
Sensitivity	0.14ng/mL
Entrez	397072
Expiration Date	Please enquire.