

## Product Datasheet

### Porcine C5a ELISA Kit (Ready to Use) (orb554331)

<b>Description</b>	Porcine Complement Component 5a (C5a) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
<b>Reactivity</b>	Porcine
<b>Range</b>	3.12-200ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Sensitivity</b>	1.12ng/mL
<b>Expiration Date</b>	Please enquire.