

Product Datasheet

Porcine CTGF ELISA Kit (Ready to Use) (orb554332)

Description	Porcine Connective Tissue Growth Factor (CTGF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection...
Reactivity	Porcine
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	F1S3P6
Sensitivity	0.242ng/mL
Entrez	100152141
Expiration Date	Please enquire.