

Product Datasheet

Porcine Cystatin C ELISA Kit (Ready to Use) (orb554335)

Description	Porcine Cystatin 3 (Cystatin C) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c...
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, urine or other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	Q0Z8R0
Sensitivity	0.056ng/mL
Entrez	733672
Expiration Date	Please enquire.