

Product Datasheet

Porcine DAO ELISA Kit (Ready to Use) (orb554338)

Description	Porcine Diamine Oxidase (DAO) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Porcine
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9TRC7
Sensitivity	0.67ng/mL
Entrez	100517436
Expiration Date	Please enquire.