

Product Datasheet

Porcine FSH ELISA Kit (Ready to Use) (orb554351)

Description	Porcine Follicle Stimulating Hormone (FSH) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Porcine
Range	6.25-400mIU/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01228
Sensitivity	2.13mIU/mL
Entrez	396895
Expiration Date	Please enquire.