

## Product Datasheet

# Porcine HMG1 ELISA Kit (Ready to Use) (orb554368)

<b>Description</b>	Porcine High Mobility Group Protein 1 (HMG1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection r...
<b>Reactivity</b>	Porcine
<b>Range</b>	62.5-4000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P12682</a>
<b>Sensitivity</b>	28.3pg/mL
<b>Entrez</b>	<a href="#">445521</a>
<b>Expiration Date</b>	Please enquire.