

Product Datasheet

Porcine IFN gamma ELISA Kit (Ready to Use) (orb554383)

Description	Porcine Interferon Gamma (IFN gamma) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Porcine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P17803
Sensitivity	5.9pg/mL
Entrez	396991
Expiration Date	Please enquire.