

Product Datasheet

Porcine IL1b ELISA Kit (Ready to Use) (orb554385)

Description	Porcine Interleukin 1 Beta (IL1b) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for...
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q29082
Sensitivity	26.2pg/mL
Entrez	396565
Expiration Date	Please enquire.