

## Product Datasheet

### Porcine IL5 ELISA Kit (Ready to Use) (orb554395)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Porcine Interleukin 5 (IL5) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...   |
| <b>Reactivity</b>        | Porcine   |
| <b>Range</b>             | 15.6-1000pg/mL  |
| <b>Form/Appearance</b>   | Ready to Use  |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.   |
| <b>Uniprot ID</b>        | <a href="#">Q9MYM5</a>  |
| <b>Sensitivity</b>       | 6.5pg/mL  |
| <b>Entrez</b>            | <a href="#">397409</a>  |
| <b>Expiration Date</b>   | Please enquire.   |