

Product Datasheet

Porcine LTF ELISA Kit (Ready to Use) (orb554398)

Description	Porcine Lactoferrin (LTF) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Porcine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P14632
Sensitivity	13.2pg/mL
Entrez	397649
Expiration Date	Please enquire.