

## Product Datasheet

### Porcine MMP9 ELISA Kit (Ready to Use) (orb554410)

<b>Description</b>	Porcine Matrix Metalloproteinase 9 (MMP9) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
<b>Reactivity</b>	Porcine
<b>Range</b>	1.56-100ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q2VI03</a>
<b>Sensitivity</b>	0.56ng/mL
<b>Entrez</b>	<a href="#">654325</a>
<b>Expiration Date</b>	Please enquire.