

Product Datasheet

Porcine MCP2 ELISA Kit (Ready to Use) (orb554412)

| | |
|--------------------------|--|
| Description | Porcine Monocyte Chemotactic Protein 2 (MCP2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ... |
| Reactivity | Porcine |
| Range | 78.1-5000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | C7F842 |
| Sensitivity | 32pg/mL |
| Entrez | 100302703 |
| Expiration Date | Please enquire. |