

Product Datasheet

Porcine NT-ProBNP ELISA Kit (Ready to Use) (orb554415)

Description	Porcine N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) Sandwich ELISA KitReady to Use: Pre-dil...
Reactivity	Porcine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	6.93pg/mL
Expiration Date	Please enquire.