

Product Datasheet

Porcine NOS3 ELISA Kit (Ready to Use) (orb554418)

Description	Porcine Nitric Oxide Synthase 3, Endothelial (NOS3) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Porcine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q28969
Sensitivity	5.4pg/mL
Entrez	397557
Expiration Date	Please enquire.