

Product Datasheet

Porcine RANTES ELISA Kit (Ready to Use) (orb554433)

Description	Porcine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) Sandwich ELISA KitR...
Reactivity	Porcine
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q2EN89
Sensitivity	1.15ng/mL
Entrez	396613
Expiration Date	Please enquire.