

## Product Datasheet

### Porcine RLN ELISA Kit (Ready to Use) (orb554434)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Porcine Relaxin (RLN) Sandwich ELISA Kit<br>Ready to Use: Pre-diluted detection reagents for convenient ... |
| <b>Reactivity</b>        | Porcine   |
| <b>Range</b>             | 31.2-2000pg/mL  |
| <b>Form/Appearance</b>   | Ready to Use  |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.             |
| <b>Sample Types</b>      | serum, plasma and other biological fluids.  |
| <b>Uniprot ID</b>        | <b>P01348</b>   |
| <b>Sensitivity</b>       | 8.36pg/mL   |
| <b>Entrez</b>            | <b>396891</b>   |
| <b>Expiration Date</b>   | Please enquire.   |