

Product Datasheet

Porcine VWF ELISA Kit (Ready to Use) (orb554464)

Description	Porcine Von Willebrand Factor (VWF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	Porcine
Range	15.6-1000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability plasma or other biological fluids.
Sample Types	plasma and other biological fluids.
Uniprot ID	Q28833
Sensitivity	4.51ng/mL
Entrez	399543
Expiration Date	Please enquire.