

## Product Datasheet

# Rabbit Apolipoprotein AI ELISA Kit (Ready to Use) (orb554470)

<b>Description</b>	Rabbit Apolipoprotein A1 (Apolipoprotein AI) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
<b>Reactivity</b>	Rabbit
<b>Range</b>	31.2-2000ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at - 20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<b>P09809</b>
<b>Sensitivity</b>	9.35ng/mL
<b>Entrez</b>	<b>100009253</b>
<b>Expiration Date</b>	Please enquire.