

Product Datasheet

Rabbit CTGF ELISA Kit (Ready to Use) (orb554483)

Description	Rabbit Connective Tissue Growth Factor (CTGF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ...
Reactivity	Rabbit
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q2WG38
Sensitivity	18.54pg/mL
Entrez	100008638
Expiration Date	Please enquire.