

Product Datasheet

Rabbit Endothelin 1 ELISA Kit (Ready to Use) (orb554486)

Description	Rabbit Endothelin 1 (Endothelin 1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents fo...
Reactivity	Rabbit
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P29560
Sensitivity	5.63pg/mL
Entrez	100009270
Expiration Date	Please enquire.