

Product Datasheet

Rabbit NT-ProBNP ELISA Kit (Ready to Use) (orb554514)

Description	Rabbit N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) Sandwich ELISA KitReady to Use: Pre-dilu...
Reactivity	Rabbit
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	11.3pg/mL
Expiration Date	Please enquire.